INSTRUCTION SHEET | ENDOSKELETAL COMPONENTS

4-Ear Socket Adapters





FIGURE 1

INSTRUCTION SHEET | ENDOSKELETAL COMPONENTS

4-Ear Socket Adapters

PRODUCT DESCRIPTION

This endoskeletal component consists of a 4-ear socket adapter. It includes a lamination cap.

KEY COMPONENTS *(Figure 1)*

4-Ear Socket Adapter

A. Titanium

B. Stainless Steel

C. Lamination Cap

This diagram is to help familiarize you with the unique parts of a 4-ear socket adapter, an endoskeletal component. These parts are referenced in the instructions and used when speaking with a technical service representative.

INTENDED USE

Endoskeletal Components are intended for aligning and attaching feet, ankles, and knees on lower limb prosthesis users.

INDICATIONS:

- Lower limb amputations

CONTRAINDICATIONS:

- None known

WARNING

- Torque evenly to specified settings.
- Do not expose this product to corrosive materials.

Failure to follow these guidelines may result in one or all of the following: noise, movement, damage, and catastrophic failure. For further questions, please contact College Park Technical Service.

WARRANTY INSPECTION/ MAINTENANCE INFORMATION

College Park recommends that you schedule your patients for check-ups according to the Warranty Inspection schedule below.

High patient weight and/or impact level may require more frequent inspections. Prosthetic components should be inspected after the first 30 days of use. We recommend you visually inspect the entire prosthesis for excessive wear and fatigue at each consultation.

WARRANTY INSPECTION SCHEDULE FOR ENDOSKELETAL COMPONENTS: 30 DAYS, THEN ANNUALLY.

TECHNICAL SPECIFICATIONS

	4-Ear Socket Adapter, w/ Pin-Lock Hole, Ti	4-Ear Socket Adapter, w/ Pin-Lock Hole, SS
Part No.	4ESAH T	4ESAH S
Material	Titanium	Titanium
Build Height	5.7 mm (0.22 in)	5.7 mm (0.22 in)
Assembly Weight	55 g	97 g
Patient Weight Limit	136 kg (300 lbs)	136 kg (300 lbs)
Warranty	3 Year, limited	3 Year, limited
Warranty	3 Year, limited	3 Year, limited
Torque (Set Screws)	15 N·m (11 ft-lbs)	15 N·m (11 ft-lbs)
Torque (Clamp Screw)	10 N·m (7.4 ft-lbs)	6 N·m (4.4 ft-lbs)

TECHNICAL ASSISTANCE / EMERGENCY SERVICE 24-7-365

College Park's regular office hours are Monday through Friday, 8:30 am – 5:30 pm (EST). After hours, an emergency Technical Service number is available to contact a College Park representative.

LIABILITY

The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer

CAUTION

College Park products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards are achieved only when College Park products are used with other recommended College Park components. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients.

CAUTION

If any problems occur with the use of this product, immediately contact your medical professional. The prosthetist and/or patient should report any serious incident* that has occurred in relation to the device to College Park Industries, Inc. and the competent authority of the Member State in which the prosthetist and/or patient is established.

*“Serious incident” is defined as any incident that directly or indirectly led, may have led, or might lead to any of the following; (a) the death of a patient, user, or other person, (b) the temporary or permanent serious deterioration of a patient's, user's, or other person's state of health, (c) a serious public health threat.

COMPLIANCE

This device has been tested according to ISO 10328 standard to three million load cycles.

Depending on patient activity this may correspond to 3-5 years of use.

ISO 10328 - LABEL

PART NUMBER	WEIGHT LIMIT (KG)	LABEL TEXT
4ESAH T	136	ISO 10328-P7-136 kg
4ESAH S		

ISO 10328 - “P” - “m”kg*)

***) Body mass limit not to be exceeded!**

For specific conditions and limitations of use, see intended use section of manufacturer's written instructions.

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©2023 College Park Industries, Inc. All rights reserved.