# INSTRUCTION SHEET | ENDOSKELETAL COMPONENTS

## Pediatric 4-Hole Pyramid Adapters and Receivers



FIGURE 1

# INSTRUCTION SHEET | ENDOSKELETAL COMPONENTS

## Pediatric 4-Hole Pyramid Adapters and Receivers

### PRODUCT DESCRIPTION

This endoskeletal component consists of a 4-hole pyramid adapter or receiver (non-rotatable or rotatable). Each includes four mounting fasteners. Pyramid receivers include four set screws. Rotatable adapters and receivers are constructed with a rotatable base, rotatable pyramid or receiver, and retainer plate.

### KEY COMPONENTS *(Figure 1)*

#### Non-Rotatable Adapters and Receivers

- A. 4-Hole Pyramid Adapter or Receiver
- B. (4) M5 x 14mm Mounting Fasteners
- C. (4) M6 x 12mm Set Screws

#### Rotatable Adapters and Receivers

- D. Rotatable Base
- E. Rotatable Pyramid or Receiver
- F. Retainer Plate
- B. (4) M5 x 14mm Mounting Fasteners
- C. (4) M6 x 12mm Set Screws

*This diagram is to help familiarize you with the unique parts of a pediatric 4-hole adapter, an endoskeletal component. These parts are referenced in the instructions and used when speaking with a technical service representative.*

### INTENDED USE

Endoskeletal Components are intended for aligning and attaching feet, ankles, and knees on lower limb prosthesis users.

#### INDICATIONS:

- Lower limb amputations

#### CONTRAINDICATIONS:

- None known

#### WARNING

- Apply Loctite® 242® on all mounting fastener threads.
- M6 x 12mm set screws do not require Loctite.
- Torque evenly to specified settings.
- Do not expose this product to corrosive materials.

*Failure to follow these guidelines may result in one or all of the following: noise, movement, damage, and catastrophic failure. For further questions, please contact College Park Technical Service.*

## WARRANTY INSPECTION/ MAINTENANCE INFORMATION

College Park recommends that you schedule your patients for check-ups according to the Warranty Inspection schedule below.

*High patient weight and/or impact level may require more frequent inspections. Prosthetic components should be inspected after the first 30 days of use. We recommend you visually inspect the entire prosthesis for excessive wear and fatigue at each consultation.*

### WARRANTY INSPECTION SCHEDULE FOR ENDOSKELETAL COMPONENTS: 30 DAYS, THEN ANNUALLY.

#### TECHNICAL SPECIFICATIONS

	4-Hole Pyramid Adapter, Standard, Ti, Pediatric	4-Hole Pyramid Adapter, Rotatable, Ti, Pediatric	4-Hole Pyramid Receiver, Ti, Pediatric	4-Hole Pyramid Receiver, Rotatable, Ti, Pediatric
<b>Part No.</b>	<b>4HPA TP</b>	<b>4HPAR TP</b>	<b>4HPR TP</b>	<b>4HPRR TP</b>
<b>Material</b>	Titanium	Titanium	Titanium	Titanium
<b>Build Height</b>	8 mm (0.31 in)	9 mm (0.35 in)	15 mm (0.59 in)	16 mm (0.63 in)
<b>Assembly Weight</b>	36 g	56 g	44 g	50 g
<b>Patient Weight Limit</b>	75 kg (165 lbs)	75 kg (165 lbs)	75 kg (165 lbs)	75 kg (165 lbs)
<b>Warranty</b>	3 Year, limited	3 Year, limited	3 Year, limited	3 Year, limited
<b>Torque (Fasteners)</b>	8 N·m (5.9 ft-lbs)	8 N·m (5.9 ft-lbs)	8 N·m (5.9 ft-lbs)	8 N·m (5.9 ft-lbs)
<b>Torque (Set Screws)</b>			8 N·m (5.9 ft-lbs)	8 N·m (5.9 ft-lbs)

## GUIDELINES

M6 x 12mm set screws with nylon patch are included with the pyramid receivers. Loctite® is not required.

M5 x 14mm fasteners included for use with all pediatric 4-hole pyramid adapters and receivers.

FASTENER SELECTION	POTENTIAL RESULT
Thread length too short	Inadequate thread engagement
Thread length too long	Bottoming-out or binding

## TECHNICAL ASSISTANCE / EMERGENCY SERVICE 24-7-365

College Park's regular office hours are Monday through Friday, 8:30 am - 5:30 pm (EST). After hours, an emergency Technical Service number is available to contact a College Park representative.

## LIABILITY

The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer

## CAUTION

College Park products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards are achieved only when College Park products are used with other recommended College Park components. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients.

## CAUTION

If any problems occur with the use of this product, immediately contact your medical professional. The prosthetist and/or patient should report any serious incident\* that has occurred in relation to the device to College Park Industries, Inc. and the competent authority of the Member State in which the prosthetist and/or patient is established.

\*'Serious incident' is defined as any incident that directly or indirectly led, may have led, or might lead to any of the following: (a) the death of a patient, user, or other person, (b) the temporary or permanent serious deterioration of a patient's, user's, or other person's state of health, (c) a serious public health threat.

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©2021 College Park Industries, Inc. All rights reserved.