# INSTRUCTION SHEET | ENDOSKELETAL COMPONENTS

## Lamination Plates




FIGURE 1

# INSTRUCTION SHEET | ENDOSKELETAL COMPONENTS

## Lamination Plates

### PRODUCT DESCRIPTION

This endoskeletal component consists of a lamination plate. It includes four dummy fasteners.

### KEY COMPONENTS *(Figure 1)*

#### Non-Threaded Plates

- |                                       |                                |
|---------------------------------------|--------------------------------|
| A. Square                             | B. Square, Pediatric           |
| C. Round, Tie-Groove                  | D. Round, Heavy Duty, Titanium |
| E. Round, Heavy Duty, Stainless Steel |                                |

#### Threaded Plates

- |                                |                    |
|--------------------------------|--------------------|
| F. Square, Threaded            | G. Round, Threaded |
| H. Round, Tie-Groove, Threaded |                    |

#### Dummy Fasteners

- |                  |                 |
|------------------|-----------------|
| I. (4) M5 x 6mm  | J. (4) M6 x 6mm |
| K. (4) M6 x 10mm |                 |

#### Recommended For Threaded Lamination Plates

- |   |
|---|
| L. Lamination Plug for 3/4"-16 Threaded Hole (not included) |
|---|

*This diagram is to help familiarize you with the unique parts of a lamination plate, an endoskeletal component. These parts are referenced in the instructions and used when speaking with a technical service representative.*

### INTENDED USE

Endoskeletal Components are intended for aligning and attaching feet, ankles, and knees on lower limb prosthesis users.

#### INDICATIONS:

- Lower limb amputations

#### CONTRAINDICATIONS:

- None known

#### WARNING

- Do not expose this product to corrosive materials.

*Failure to follow these guidelines may result in one or all of the following: noise, movement, damage, and catastrophic failure. For further questions, please contact College Park Technical Service.*

## INSTRUCTIONS

1. Apply resin deterrent to lamination fastener threads and finger-tighten to plate. Then, apply resin deterrent to cap screw openings.
2. Attach plate to cast or first-stage laminate using method of choice (i.e. adhesive, lamination plug, lock dummy, etc).
3. Continue with fabrication method of choice.
4. Once cured, grind distal end of socket flush to heads of lamination fasteners, then remove fasteners and discard.
5. Attach 4-hole adapter of choice. Apply Loctite® 242® to mounting fastener threads and torque to 12 N·m (8.8 ft-lbs.)]

## TECHNICAL ASSISTANCE / EMERGENCY SERVICE 24-7-365

College Park's regular office hours are Monday through Friday, 8:30 am – 5:30 pm (EST). After hours, an emergency Technical Service number is available to contact a College Park representative.

## LIABILITY

The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer

## CAUTION

College Park products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards are achieved only when College Park products are used with other recommended College Park components. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients.

## CAUTION

If any problems occur with the use of this product, immediately contact your medical professional. The prosthetist and/or patient should report any serious incident\* that has occurred in relation to the device to College Park Industries, Inc. and the competent authority of the Member State in which the prosthetist and/or patient is established.

\*'Serious incident' is defined as any incident that directly or indirectly led, may have led, or might lead to any of the following: (a) the death of a patient, user, or other person, (b) the temporary or permanent serious deterioration of a patient's, user's, or other person's state of health, (c) a serious public health threat.

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©2021 College Park Industries, Inc. All rights reserved.