

**college park**  
TECHNOLOGY for the HUMAN RACE

## technical instructions

التعليمات الفنية • Teknisk vejledning • Technische Anleitung • Τεχνικές οδηγίες • Instrucciones técnicas  
Tekniset ohjeet • Instructions techniques • הוראות טכניות • Istruzioni tecniche • Technische instructies  
Teknische instruksjoner • Instrukcje techniczne • Instruções técnicas • Instruções técnicas  
Технические инструкции • Technické pokyny • Tekniska anvisningar • Teknik Talimatlar • 技术说明

**PACKAGE CONTENTS**

(1) College Park Aura Knee

(1) Knee External Lock Kit

This diagram is to help familiarize you with the unique parts of the Aura Knee. These parts are referenced in the instructions and used when speaking with a technical service representative.

**KEY COMPONENTS (FIGURE 1)****A.** Pyramid**B.** Upper Assembly**C.** Lower Assembly**D.** 30mm Pylon Receiver**E.** 4mm Clamp Screw (Torque 10 N-m)**F.** Release Lever**G.** Arm Linkage**H.** Lock Button**I.** External Lock Kit**J.** Knee Cap**K.** Adjustment Cover**PRODUCT DESCRIPTION**

The Aura Knee is constructed with an integrated pyramid (proximal) and 30 mm pylon receiver (distal).

**INTENDED USE**

The Aura Knee, intended for transfemoral amputees, is a prosthetic device designed to restore some function of an anatomical knee joint.

**⚠ INDICATIONS:**

Lower limb amputations above the knee

**⚠ CONTRAINDICATIONS:**

None known

**TECHNICAL SPECIFICATIONS****Material** Aluminum, Stainless Steel**Assembly Weight** 665 g**Joint Type** Polycentric**Patient Weight Limit** 300 lb (136 kg)**Flexion** 135°**Torque (Clamp Screw)** 10 N·m**Build Height** 0.65 in (1.7 cm)**TOOLS RECOMMENDED**

(1) 4 mm Hex Key

**BUILD HEIGHT (FIGURE 2)****A** Overall height **6.32 in (16.1 cm)****B** Dome to knee center **0.65 in (1.7 cm)****C** Dome to tube end contact **4.32 in (11.0 cm)****BENCH ALIGNMENT (FIGURE 3)****A.** Determine the heel height**B.** Determine the socket flexion**C.** The load line bisects through the pylon**STATIC ALIGNMENT**

With the alignment reference line through the pylon, plantarflex or dorsiflex the foot until the load line is balanced between 1/3 heel and 2/3 toe lever.

**KNEE ALIGNMENT****More Stable** = slide the knee posterior**More Dynamic** = slide the knee anterior**DYNAMIC ADJUSTMENTS**

All adjustments can be made using a 4mm Allen wrench.

Note: Patient must be sitting when adjustments are made.

Note: The effect of a dynamic adjustment can be felt with as little as 1/8 - 1/4 turn.

Caution: Flexion and extension must be possible at all settings.

## STANCE PHASE

### STABILITY ADJUSTMENT (FIGURE 4)

**Factory setting = Screw is at minimum (most stable). Should only be adjusted a maximum of 2 turns.**

Turning clockwise straightens the knee, causing it to be closer to the point where it will trigger or flex freely when weight is applied.

Note: This adjustment can be used to remove play to the lock or to require weight/toe load to lock. Adjustment may cause the knee to be difficult or unable to lock.

SYMPTOM	DESIRED RESULT	SCREW ADJUSTMENT
Alignment too stable; Swing too difficult to initiate	Decrease Knee Stability	Turn screw clockwise
Alignment too unstable; Swing too easy to initiate	Increase Knee Stability	Turn screw counterclockwise

### STANCE FLEXION ADJUSTMENT (FIGURE 5)

**Factory setting = Screw is at minimum (no stance flexion).**

Controls the amount of stance flexion. Adjustments can be made to increase the amount of stance flexion up to 6 degrees

SYMPTOM	DESIRED RESULT	SCREW ADJUSTMENT
Not enough flexion	Increase Stance Flexion	Turn screw counterclockwise
Too much flexion	Decrease Stance Flexion	Turn screw clockwise

## SWING PHASE

### FRICTION SWING (FIGURE 6)

**Factory setting = Screw is at minimum.**

Controls the rate of heel rise.

SYMPTOM	DESIRED RESULT	SCREW ADJUSTMENT
Excessive heel rise	Increase Friction	Turn screw clockwise
Heel rise not sufficient	Decrease Friction	Turn screw counterclockwise

### EXTENSION ASSIST (FIGURE 7)

**Factory setting = Screw is at minimum.**

Controls the rate of extension assist.

SYMPTOM	DESIRED RESULT	SCREW ADJUSTMENT
Extension too slow	Increase Extension Assist	Turn screw clockwise
Extension too fast	Decrease Extension Assist	Turn screw counterclockwise

## LOCKING OPTIONS

**Factory setting = Normal Mode**

The College Park Aura can be converted to a locking knee joint, allowing the patient to control when the knee is locked or unlocked.

<b>Normal Mode</b>	Lock is disabled; Knee utilizes polycentric geometry for stability during stance.
<b>Locking Mode</b>	Lock is enabled; Allows for manual locking and unlocking of knee

## LOCKING THE KNEE (FIGURE 8)

1. Push the lock button from the left side.

Caution: Test the knee to ensure it does not flex, prior to patient fitting.

2. Continue with installation of the External Lock Kit (optional).

## UNLOCKING THE KNEE

Push the lock button from the right side. The knee will now be in normal mode..

Caution: Test the knee to ensure it swings freely without weight applied, prior to patient fitting.

## USING THE EXTERNAL LOCK KIT (FIGURE 9)

To operate, the knee must be in locking mode, and the external lock kit must be installed. Refer to External Lock Kit Fabrication Instructions for detailed information and mounting instructions.

Caution: Test the dual function lever after installation and all adjustments to ensure proper function.

Note: Instruct the patient on proper use of the dual function lever. Make sure they understand the operating positions. Visit [www.college-park.com](http://www.college-park.com) to find fabrication instructions.

<b>A. LOCKED POSITION</b>	Standing/Walking: The knee will not flex.
<b>B. LOCK RELEASE POSITION</b>	Pull upward on the lever only half-way, without passing the detents. Releasing the lever will return it to the locked position.
<b>C. UNLOCKED POSITION</b>	Pull upward on the lever all the way, passing the detents. There will be a “click” sound and feeling. Push down on the lever to return to the locked position.

## TO SIT DOWN

Patient should pull upward on the dual function lever to allow the knee to flex.

## USE IN WATER

- The Aura Knee has been approved for use in water.
- The knee should be switched into lock mode when using it in or near water.
- Use caution when walking on wet surfaces.
- Always dry knee thoroughly after exposure to moisture. Rinse and dry if exposed to salt water or chlorinated water.

## ⚠️ WARNINGS

- Flexion and extension must be possible at all settings.
- The knee is not intended for kneeling or squatting.
- Avoid pinching hazards! Do not place fingers near the flexing area of the knee.
- Patient concerns about the function should be reported to the prosthetist immediately, including but not limited to noise, sudden loss of function, brake release sticking, etc.
- Do not disassemble the knee. Contact College Park to arrange a repair or replacement.
- Contaminants such as dirt and the use of lubricants or powder may affect the function of the knee brake and lead to failure.
- Do not use compressed air to clean the knee as it can push dirt inside of the knee.
- If the locking function is utilized, verify that the manual lever and nylon lanyard are not obstructed by the cosmesis.
- If modifying the prosthesis, ensure that the nylon lanyard is secure to prevent damage from grinding tools.

Failure to follow these technical instructions or use of this product outside the scope of its Limited Warranty may result in injury to the patient or damage to the product.

## WARRANTY INSPECTION AND MAINTENANCE INFORMATION

College Park recommends that you schedule your patients for check-ups per the warranty inspection schedule below.

High patient weight or activity level may require more frequent inspections. We recommend you visually inspect the following applicable parts for excessive wear and fatigue at each warranty inspection.

- Knee Assembly, External Lock

**WARRANTY INSPECTION SCHEDULE FOR COLLEGE PARK AURA: SIX MONTHS, THEN ANNUALLY.**

## TECHNICAL ASSISTANCE / EMERGENCY SERVICE 24-7-365

College Park's regular office hours are Monday through Friday, 8:30 am – 5:30 pm (EST). After hours, an emergency Technical Service number is available to contact a College Park representative.

### LIABILITY

The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer

#### CAUTION

If the product is exposed to conditions outside of the recommended environmental conditions inspect for signs of changes or loss of functionality during use which may include poor response, lack of swing flexion initiation, abnormal noises, etc. If there are noticeable changes or loss of function discontinue use of the product.

#### CAUTION

College Park products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards are achieved only when College Park products are used with other recommended College Park components. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients.

#### CAUTION

If any problems occur with the use of this product, immediately contact your medical professional. The prosthetist and/or patient should report any serious incident\* that has occurred in relation to the device to College Park Industries, Inc. and the competent authority of the Member State in which the prosthetist and/or patient is established.

\*“Serious incident” is defined as any incident that directly or indirectly led, may have led, or might lead to any of the following; (a) the death of a patient, user, or other person, (b) the temporary or permanent serious deterioration of a patient’s, user’s, or other person’s state of health, (c) a serious public health threat.

### COMPLIANCE

This device has been tested according to ISO 10328 standard to three million load cycles.

Depending on patient activity this may correspond to 3-5 years of use.

#### ISO 10328 - LABEL

##### WEIGHT LIMIT (KG)

136

##### LABEL TEXT

ISO 10328-P6-136 kg

ISO 10328 - “P” - “m”kg\*) 

\*) *Body mass limit not to be exceeded!*  
For specific conditions and limitations of use, see intended use section of manufacturer's written instructions.

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©2024 College Park Industries, Inc. All rights reserved.