



## WARRANTY INSPECTION AND MAINTENANCE INFORMATION

College Park recommends that you schedule your patients for check-ups according to the warranty inspection schedule below.

High patient weight and/or impact level may require more frequent inspections. We recommend you visually inspect the following applicable parts for excessive wear and fatigue at each warranty inspection.

- Composite and Adapter
- Foot Shell

WARRANTY INSPECTION SCHEDULE FOR TEMPO : ANNUALLY.

## TECHNICAL ASSISTANCE / EMERGENCY SERVICE 24-7-365

College Park's regular office hours are Monday through Friday, 8:30 AM – 5:30 PM (EST). After hours, an emergency Technical Service number is available to contact a College Park representative.



**COLLEGE PARK INDUSTRIES, INC**  
27955 College Park Dr.  
Warren, MI 48088 USA

**EMERGO EUROPE**  
Prinsessegracht 20, 2514 AP  
The Hague, Netherlands


**CE**  
**MADE IN THE USA**  
311 INS TO TIS 171208

**Australian Sponsor**  
**EMERGO AUSTRALIA**  
Level 20, Tower II  
Darling Park  
201 Sussex Street  
Sydney, NSW 2000  
Australia

©2018 College Park Industries, Inc. All rights reserved. College Park, Tempo, and Gait Matching are registered trademarks of College Park Industries. Intelliweave is a trademark of College Park Industries.



*technical  
instructions*



This diagram is to help familiarize you with the unique parts of the Tempo. These parts are referenced in the instructions and used when speaking with a technical service representative.

## KEY COMPONENTS

1. Integrated Pyramid
2. Housing (Ankle)
3. Bonded Foot Shell

## PACKAGE CONTENTS

(1) Tempo foot, with bonded foot shell

## TOOLS RECOMMENDED

4mm Hex Key

## WARNING

- Do not expose this product to pH extremes or corrosive materials (water, salt water or other liquids).
- Failure to follow these technical instructions or use of this product outside the scope of its Limited Warranty may result in injury to the patient or damage to the product.
- Disassembly, modification or removal from the foot shell, or any further disassembly or modification of components will void the warranty.

## PROTECTIVE COVER ON DOME

Remove the protective cover on dome after alignment is completed and before patient leaves clinic.

## MOUNTING

Use only high quality proximal endoskeletal components (30 mm). The service-free Trés foot comes permanently assembled in the foot shell.

## STATIC ALIGNMENT

For optimal function have the patient's weight balanced between the forefoot and heel. The load line divides the foot at the 1/3 heel lever and 2/3 toe lever.

## DYNAMIC ADJUSTMENTS

DESIRED RESULT	ALIGNMENT CHANGE
FIRMER TOE RESPONSE	Plantarflex the Tempo or move load line posterior
SOFTER TOE RESPONSE	Dorsiflex the Tempo or move load line anterior
FIRMER HEEL RESPONSE	Dorsiflex the Tempo or move load line anterior
SOFTER HEEL RESPONSE	Plantarflex the Tempo or move load line posterior

## GAIT MATCHING® GUIDELINES

Contact College Park Technical Service if you have questions about modifying the gait match.