college park
TECHNOLOGY *for the HUMAN RACE*

technical instructions

التعليمات الفنية · Teknisk vejledning · Technische Anleitung
Τεχνικές οδηγίες · Instrucciones técnicas · Tekniset ohjeet
Instructions techniques · הוראות טכניות · Istruzioni tecniche
Technische instructies · Tekniske instruksjoner · Instrukcje techniczne
Instruções técnicas · Instruções técnicas · Технические инструкции ·
Technické pokyny · Tekniska anvisningar · Teknik Talimatlar · 技术说明

FIGURE 1



FIGURE 2

Package Contents

(1) Trustep Foot	(1) Foot Shell
(1) CPI Sock	(1) Sealing Boot (upon request)
(1) TruLube Lubricant	

Tools Recommended

(1) #2 Phillips Screwdriver	(1) Torque Wrench
(1) 4 mm and 6 mm Hex Keys	(1) Foot Horn
(1) Pin Guide	(1) Preloader Tool

This diagram is to help familiarize you with the unique parts of the Trustep. These parts are referenced in the instructions and used when speaking with a technical service representative.

Key Components (Figure 1)

A. Ankle Bone	B. Heel Bone	C. Forefoot Bone	D. Rear Bumper
E. Midstance Pad	F. Front Bumper	G. Ankle Bushings (2)	H. Lower Bushings (2)
I. Axial Pins (2)	J. Axial Pin Screws (2) Torque 4 N-m (36 in-lbs)	K. Preload Insert Torque 2 N-m (18 in-lb)	L. Insert Screw

• CPI Sock (not shown) • Foot Shell (not shown)

Optional EXO Mount

M. Exo Pyramid Tool (additional purchase)	N. Exo Block Kit (additional purchase)	O. Exo Ankle Bone	P. Exo Mounting Bolt (included w/exo block kit) Torque 61 N-m (45 ft-lbs)
---	--	--------------------------	--

PRODUCT DESCRIPTION:

This prosthetic foot device is constructed with an endo or exo ankle bone option, heel bone, forefoot bone, bumpers, bushings, axial pins, preload insert, and fasteners. The ankle, heel, and forefoot bones are assembled with the axial pins and screws. The forefoot is secured with the preload insert and screw.

INTENDED USE:

The Trustep is a prosthetic foot designed to replace one or more functions of the biologic human foot.

⚠ INDICATIONS:

- Lower limb amputations

⚠ CONTRAINDICATIONS:

- None known

⚠ PROTECTIVE COVER ON DOME

Remove the protective cover on dome after alignment is completed and before patient leaves clinic.

TECHNICAL SPECIFICATIONS

FOOT SIZE	MOUNTING	WEIGHT LIMIT	BUILD HEIGHT	FOOT WEIGHT*
22-31 cm	Endo	330 lbs / 150 kg	3.5 in / 8.7 cm	678 g
	Exo	352 lbs / 160 kg	6.0 in / 15.2 cm	

*26 cm foot with shell

GAIT MATCHING® GUIDELINES

The gait match determines the firmness of the foot based on the user's specifications (foot size, patient weight, and impact level).

FIRMNESS CATEGORIES

Refer to the chart below to determine the correct firmness category.

Note: Incorrect category selection may result in poor device function. Contact College Park Technical Service if you have questions about category selection.

WEIGHT LBS	0-74	75-100	101-120	121-140	141-160	161-180	181-200	201-220	221-250	251-275	276-300	301-330	331-352
WEIGHT KG	0-33	34-45	46-54	55-63	64-72	73-81	82-91	92-100	101-113	114-125	126-136	137-150	151-160
SIZE CM	22-31												EXO ONLY
LOW IMPACT	1	2	3	4	4	5	6	6	7	7	8	8	9
MODERATE IMPACT	1	2	3	4	5	6	6	7	7	7	8	9	9
HIGH IMPACT	2	3	4	5	6	6	7	7	8	8	8	9	9

GAIT MATCH CHART

Components may be changed as needed. Determine the user's firmness category, then refer to the chart below for the recommended components.

For gait match adjustments, see section: Dynamic Adjustments.

COMPONENT	FIRMNESS CATEGORY								
	1	2	3	4	5	6	7	8	9
FRONT BUMPER	1	2	3	4	5	6	7	8	8
ANKLE BUSHINGS*	S	S	S	S	M	M	M	F	F
MIDSTANCE PAD*	M	M	M	M	M	M	M	F	F
REAR BUMPER	XS	XS	S	S	M	H	1X	2X	3X
LOWER BUSHINGS	ONE FIRMNESS	ONE FIRMNESS	ONE FIRMNESS	ONE FIRMNESS	ONE FIRMNESS	ONE FIRMNESS	ONE FIRMNESS	ONE FIRMNESS	ONE FIRMNESS

ENDOSKELETAL MOUNTING

Use only high quality proximal endoskeletal components (30 mm).

EXOSKELETAL MOUNTING

1. Remove Exo Ankle Bone, then attach to College Park Exo Block with anti-rotation pin(s) oriented. Apply Loctite® 242 to mounting bolt. Torque to 61 N·m (45 ft·lbs). To skip alignment and lamination, go to Step 9.
2. If using the align-able Exo option, attach Exo Pyramid Tool to Exo Block with four 6 mm screws. Torque to 15 N·m (11 ft·lbs).
3. Attach 30 mm endo components to the Exo Pyramid Tool and temporarily mount the socket.
4. Re-attach the ankle bone to foot, donning CPI Sock and foot shell, then perform a dynamic alignment.
5. Remove foot from Exo Block.
6. Mount aligned prosthesis in transfer jig. Lock socket and Exo Block in position.
7. Remove endo components and Exo Pyramid Tool.
8. Use desired method to span Exo Block to the socket and remove from jig. Shape and laminate to desired finish. Do not remove foam from the top of the Exo Block
9. Re-attach ankle bone to heel bone. Re-assemble foot, don CPI Sock and foot shell.

See the Trustep Assembly / Disassembly videos on our YouTube channel!

ANKLE REMOVAL

(ACCESS TO REAR BUMPER, ANKLE BUSHINGS AND MIDSTANCE PAD)

1. Use the Foot Horn to don/doff the foot shell. Remove the CPI Sock and replace as needed.
2. Use a 6 mm and 4 mm hex keys to remove the axial pin screw.
3. Attach the pin guide to the axial pin. Push and rotate clockwise.
4. Unscrew the axial pin from the pin guide and remove it from the ankle bone.
5. Remove the ankle bone to access the rear bumpers, ankle bushings and midstance pad.
6. For re-assembly, lubricate the axial pin, inside and outside of the ankle bushings, ankle bushing pocket and the contact surface of the midstance pad, then reverse steps 1-4. Torque the axial pin screw to 4 N·m (36 in·lbs).

DO NOT lubricate front and rear bumpers.

FOREFOOT REMOVAL

(ACCESS TO FRONT BUMPER AND LOWER BUSHINGS)

1. To remove the Preload Insert, position the Preloader Tool on the foot to remove the load off of the insert. Tighten the handle down on the forefoot then swing the cam lever to either side.
For detailed instructions refer to the Preloader Tool Instruction Sheet.
2. Remove the front bumper and axial pin as instructed (above) for the Ankle Bone.
3. Remove the forefoot bone to access the lower bushings.
4. For re-assembly, lubricate the axial pin, inside and outside of the lower bushings, then reverse steps 1-4.
Torque the axial pin screw to 4 N·m (36 in·lbs)

STATIC ALIGNMENT (Figure 2)

For optimal function, balance the patient's weight evenly between the heel and toe.

- The Trustep was designed with a 1/2" (13 mm) heel rise.
- The load line divides the foot at 1/3 heel lever and 2/3 toe lever.

DYNAMIC ADJUSTMENTS

DESIRED RESULT	ALIGNMENT CHANGE	COMPONENT CHANGE
FIRMER TOE RESPONSE	Plantarflex the Trustep or move load line posterior	Front Bumper up one step firmer
SOFTER TOE RESPONSE	Dorsiflex the Trustep or move load line anterior	Front Bumper down one step softer or use MIN Insert
FIRMER HEEL RESPONSE	Dorsiflex the Trustep or move load line anterior	Rear Bumper up one step firmer
SOFTER HEEL RESPONSE	Plantarflex the Trustep or move load line posterior	Rear Bumper down one step softer
INCREASE/DECREASE ROTATION, INVERSION OR EVERSION	—	Change Ankle Bushings

WARNING

- Do not expose this product to pH extremes or corrosive materials (water, salt water or other liquids).
- Failure to follow these technical instructions or use of this product outside the scope of its Limited Warranty may result in injury to the patient or damage to the product.
- Disassembly, modification or removal from the foot shell, or any further disassembly or modification of components will void the warranty.

RESIDUAL RISK STATEMENT

NOTICE OF RESIDUAL RISK

During fitting process, ensure that CPI sock does not become pinched between foot and endoskeletal componentry.

WARRANTY INSPECTION / MAINTENANCE INFORMATION

College Park recommends that you schedule your patients for check-ups according to the warranty Inspection schedule below.

High patient weight and/or impact level may require more frequent inspections. Soft component wear depends on the patient weight, impact level and environment. We recommend you inspect the following applicable parts for excessive wear and fatigue at each warranty inspection and replace as needed.

- Soft components (disassemble, inspect & re-lubricate)
- CPI Sock
- Composites and Adapters
- Foot Shell

WARRANTY INSPECTION SCHEDULE FOR TRUSTEP: SIX MONTHS, THEN ANNUALLY.

TECHNICAL ASSISTANCE / EMERGENCY SERVICE 24-7-365

College Park's regular office hours are Monday through Friday, 8:30 am – 5:30 pm (EST). After hours, an emergency Technical Service number is available to contact a College Park representative.

LIABILITY

The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer

CAUTION

College Park products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards are achieved only when College Park products are used with other recommended College Park components. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients.

CAUTION

If any problems occur with the use of this product, immediately contact your medical professional. The prosthetist and/or patient should report any serious incident* that has occurred in relation to the device to College Park Industries, Inc. and the competent authority of the Member State in which the prosthetist and/or patient is established.

**"Serious incident" is defined as any incident that directly or indirectly led, may have led, or might lead to any of the following; (a) the death of a patient, user, or other person, (b) the temporary or permanent serious deterioration of a patient's, user's, or other person's state of health, (c) a serious public health threat.*

COMPLIANCE

This device has been tested according to ISO 10328 standard to two million load cycles. Depending on patient activity this may correspond to 2-3 years of use.

ISO 10328 - LABEL

FOOT SIZE	WEIGHT LIMIT (KG)	LABEL TEXT
22-31 cm (Endo)	150	ISO 10328-P7-150 kg
22-31 cm (Exo)	166	ISO 10328-P8-166 kg

ISO 10328 - "P" - "m"kg*

***) Body mass limit not to be exceeded!**
For specific conditions and limitations of use, see intended use section of manufacturer's written instructions.

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©2023 College Park Industries, Inc. All rights reserved.