



TERRAIN iW PART NUMBER

	SIDE	SHELL COLOR		SIZE		FIRMNESS CATEGORY
TN iW	L OR R	C, B, OR T	EN	21-30	D	1-4
	(L) Left, (R) Right	(C) Caucasian, (B) Brown, (T) Tan				



NEW

College Park designs and manufactures a full line of anatomically correct, custom built prosthetic solutions. College Park is ISO 13485 certified with an extensive focus on quality testing, and proudly sources and produces all products entirely in the USA. With precise engineering and quality manufacturing, College Park is committed to innovating human locomotion solutions for users all over the world.





Patients have different goals and activity levels. It's our mission to create prosthetic solutions that maximize benefits and pinpoint common interests. With a focus on low-impact users, College Park combines some of its most popular and trusted technology in the new Terrain iW. This lightweight foot offers tri-axial movement on an Intelliweave composite base. Featuring anatomical ankle motion, flexible springs, and a durable heel wedge, the Terrain iW was designed for comfort and stability.

TRI-AXIAL ANKLE MOVEMENT

Building upon the success of the popular Terrain foot, the Terrain iW uses a similar ankle to achieve movement on all three planes—the sagittal, coronal, and transverse. Whether you're relaxing with family, grocery shopping, or simply walking to get the mail, this foot can increase ground compliance for everyday obstacles.

TRANSVERSE ROTATION



INVERSION/EVERSION



PLANTAR/DORSIFLEXION



9.1° of total motion*

15.8° of total motion*

*including movement in springs



SUPPORTIVE AND SERVICE-FREE

User friendly for patients and clinicians, the Terrain iW is service-free! With permanent bumpers inside the ankle housing, there's no need to disassemble the foot and replace urethane parts. The Terrain iW's durable heel wedge is designed to enhance stability upon heel strike—which may improve ground contact while walking.



iW: INTELLIWEAVE

Intelliweave technology aims to provide flexibility and strength. These fibers are woven by hand in a 3D pattern and gait-matched according to a patient's unique specifications. As the Terrain iW responds to a user's movement, there may be a lighter impact on the limb and socket.



TECHNICAL SPECIFICATIONS

FOOT BASE

Intelliweave®

MOUNTING

Endo

SIZES

21-30 cm

WEIGHT LIMIT

21-24 cm | 220lbs (100kg)
25-30 cm | 300lbs (136kg)

BUILD HEIGHT

2.5 in (21-24 cm)
2.7 in (25-26 cm)
2.8 in (27-30 cm)

WATER CAPABILITY: Waterproof in Fresh Water

WEIGHT*

645 g

WARRANTY

2 years

L CODES**

L5972 & L5986

PDAC APPROVED!

IMPACT LEVEL



ENVIROSHELL®†



*26 cm model, including pyramid and Enviroshell. **The listing of codes with these products should not be construed as a guarantee for coverage or payment. Ultimate responsibility for the coding of services/products rests with the individual practitioner.
† Caucasian, Tan, or Brown.

