

ICON™

microprocessor knee
user friendly design
streamlined setup process
responses in real time



The Icon microprocessor knee offers a user-friendly experience for patients and clinicians. Featuring responsive sensors, streamlined setup, and the intuitive Stride Studio app, this dynamic knee was designed to enhance user confidence and trust. Clinicians can easily access all adjustments, custom modes, and more in their version of Stride Studio. The Icon has no service requirement with College Park's remote updates and unique warranty process. Our engineers have developed patented technology to put the Icon in a league of its own—improving angular accuracy, response time, and weight savings. Equipped with a long-lasting battery and IP68 rating, the Icon is the versatile solution for moderate to high activity users.

PART NUMBER










IC MPKA 1000

RECOMMENDED ADAPTERS



	PART NUMBER	DESCRIPTION
	TA T250 30	250 mm Low-Profile Tube Adapter, 30 mm diameter, Titanium
	TA T430 30	430 mm Low-Profile Tube Adapter, 30 mm diameter, Titanium
	TC T 30	Tube Clamp, 30 mm diameter, Titanium
	3PA R	3-Prong anchor, rotatable

K2-4 ICON

REPLACEMENT PARTS

	PART NUMBER	DESCRIPTION
	IC WAR3	Icon Microprocessor Knee 3 Year Extended Warranty
	IC KPK	Icon Knee Pad Replacement Kit
	IC BATT	Icon Knee Replacement Battery
	IC CHG US	Icon Knee Charging Kit, US
	IC CHG EU	Icon Knee Charging Kit, Euro
	IC CHG UK	Icon Knee Charging Kit, UK
	IC CB US	Icon Knee Charging Blade, US
	IC CB EU	Icon Knee Charging Blade, EU
	IC CB UK	Icon Knee Charging Blade, UK
	IC CB AU	Icon Knee Charging Blade, AU

HARDWARE & TOOLS

	PART NUMBER	DESCRIPTION
	LPTC T	Low Profile Threaded Connector, Ti
	CA GELK	Knee Cyanoacrylate Adhesive Gel
	CP HK 4	Hex Key, 4mm

TECH SPECS

DISTAL ATTACHMENT

Pyramid

WEIGHT

1530 g

WEIGHT LIMIT

330 lbs | 150 kg

FLEXION RANGE

130°

OVERALL HEIGHT

10.5 in | 267 mm

DOMES TO KNEE CENTER

0.67 in | 17 mm

WARRANTY

3 years

L-CODES*

L5828, L5845, L5848,
L5850, L5856, L5925

WATER CAPABILITY



FRESH WATER PROOF

*The listing of codes with these products should not be construed as a guarantee for coverage or payment. Ultimate responsibility for the coding of services/products rests with the individual practitioner.