

# VICTOR™

## pneumatic knee

easy adjustments  
stance flexion  
variable cadence



Step victoriously. The Victor knee is a single axis, pneumatic knee joint featuring a friction brake to control stance phase. This knee is the ideal solution for low to moderate activity users, as it can adjust to accommodate both slow and fast paced walking. It also works well for patients looking to increase their activity level after rehabilitation. The Victor is a safe, economical option that instills confidence and freedom in patients.

### PART NUMBER



# VPKA 200

### HARDWARE & TOOLS

### TECH SPECS

<b>MATERIAL</b> Aluminum	<b>WEIGHT</b> 710 g
<b>WEIGHT LIMIT</b> 275 lbs	<b>FLEXION RANGE</b> 145°
<b>WARRANTY</b> 2 years	<b>L-CODES*</b> L5822, L5845, L5850

\*The listing of codes with these products should not be construed as a guarantee for coverage or payment. Ultimate responsibility for the coding of services/products rests with the individual practitioner.

	PART NUMBER	DESCRIPTION
	CP KSB	Knee Stop Bumper Replacement Kit
	CP HK 4	Hex Key, 4mm