


AURA™

polycentric knee

simple 4mm adjustments
unique manual lock
waterproof

K2



The Aura polycentric knee was designed to offer stability and relief to low-impact users. Featuring a unique push-button lock and simple 4mm adjustments, the waterproof Aura goes beyond mere physical support. It opens doors for mobility goals and newfound peace of mind.

PART NUMBER

APKA 150

TECH SPECS


DISTAL ATTACHMENT	WEIGHT
30mm tube clamp	665 g
WEIGHT LIMIT	FLEXION RANGE
300 lbs 136 kg	135°
OVERALL HEIGHT	DOME TO KNEE CENTER
6.32 in 161 mm	0.65 in 17 mm
WARRANTY	L-CODES*
2 years	L5611, L5845, L5850, L5925

WATER CAPABILITY








WATERPROOF




*The listing of codes with these products should not be construed as a guarantee for coverage or payment. Ultimate responsibility for the coding of services/products rests with the individual practitioner.



RECOMMENDED ADAPTERS

	PART NUMBER	DESCRIPTION
	TA T250 30	250 mm Low-Profile Tube Adapter, 30 mm diameter, Titanium
	TA T430 30	430 mm Low-Profile Tube Adapter, 30 mm diameter, Titanium
	TC T 30	Tube Clamp, 30 mm diameter, Titanium
	3PA R	3-Prong anchor, rotatable
	LPTCT	Low Profile Threaded Connector, Ti

HARDWARE & TOOLS

	PART NUMBER	DESCRIPTION
	CP HK 4	Hex Key, 4mm
	CP KELK	Knee External Lock Kit (included with knee)
	CP LK	Lanyard Kit (included with knee)