



AVIATOR[®]

vertical shock absorption
high-performance comfort
material fusion



Soar through life with the Aviator. This new, lightweight foot is designed to provide high energy return, comfort, and vertical shock absorption in every step. Combining flexible Intelliweave[®] composites with dynamic carbon fiber, the sleek Aviator uses a hybrid material approach in its springs. These features offer superior range of motion to accommodate a wide range of activities and users. When life gets turbulent, the Aviator can help you stay the course.

TECH SPECS

FUNCTION

Vertical Shock

BUILD HEIGHT

6.5 in

MOUNTING

Endo (IP)

WEIGHT LIMIT

21-25 cm | 275 lbs

26-30 cm | 330 lbs

L-CODE[‡]

L5987

PDAC APPROVED!

WATER CAPABILITY



FRESH WATER PROOF

SIZES

21-30 cm

WEIGHT*

979 g

IMPACT LVL



WARRANTY

3 years

ENVIROHELL[†]



standard

*26 cm with Enviroshell. ‡ The listing of codes with these products should not be construed as a guarantee for coverage or payment. Ultimate responsibility for the coding of services/products rests with the individual practitioner. †Caucasian, Tan, Brown and Jet Black.

AVIATOR PART NUMBER

	SIDE	SHELL COLOR*		SIZE		FIRMNESS CATEGORY
AV	L OR R	C, T, B OR J	EN	21-30	D	1-5

* (C)aucasian, (T)an, (B)rown or (J)et Black. ‡ (N)arrow or (W)ide. Example part number: AV LBEN 21 D1

AVIATOR FIRMNESS CATEGORIES

PATIENT WEIGHT	0-140lbs	141-180	181-220	221-275	276-330
SIZE	21-25cm				
MODERATE-HIGH IMPACT	1	2	3	4	
SIZE	26-30cm				
MODERATE-HIGH IMPACT	1	2	3	4	5

ENVIROHELL PART NUMBER

	SIDE	SHELL COLOR**	SIZE
HZ S	L OR R	C, T, B, OR J	21-30

** (C)aucasian, (T)an, (B)rown, or (J)et Black



Enviroshell included with every foot order


INCLUDED

Endo foot includes: foot shell, CPI sock and wedge kit

K3/4



AVIATOR



HARDWARE

	PART NUMBER	DESCRIPTION
	VY HWK 2125	Voyager Wedge Kit, 21-25
	VY HWK 2630	Voyager Wedge Kit, 26-30cm
	CP SSET	Super Set Screws (8pc) M8 x 14mm

TOOLS

	PART NUMBER	DESCRIPTION
	TFH 800	Foot Horn, Stainless Steel
	CP HK 4	Hex Key, 4mm
	CA 409	Cyanoacrylate Adhesive, 409 Gel (CA glue)

SHELLTREAD

	PART NUMBER (LEFT/RIGHT)	DESCRIPTION	ADDITIONAL INFORMATION
	HZ TLN xx HZ TRN xx	Shelltread (left or right) Neutral (indicate size, width)	Increases durability and prevents premature wear of foot shell.
	HZ TLJ xx HZ TRJ xx	Shelltread (left or right) Jet Black (indicate size, width)	

SOCKS

	PART NUMBER	DESCRIPTION	ADDITIONAL INFORMATION
	HZ SS 2125	CPI Sock, 21-25cm	Highly durable sock that acts as a noise barrier and protects the foot from debris.
	HZ SS 2630	CPI Sock, 26-30cm	
	HZ SSK 2125	CPI Sock, 6-pack, 21-25cm	
	HZ SSK 2630	CPI Sock, 6-pack, 26-30cm	