

TERRAIN® iW

tri-axial technology
intelliweave comfort
service-free solution



TECH SPECS

FUNCTION

Tri-axial

BUILD HEIGHT

(21-24 cm) 2.5 in
(25-26 cm) 2.7 in
(27-30 cm) 2.9 in

MOUNTING

Endo (IP)

WEIGHT LIMIT

21-24 cm | 220 lbs
25-30 cm | 300 lbs

L-CODE‡

L5972
L5986

WATER CAPABILITY



FRESH WATER PROOF

SIZES

21-30 cm

WEIGHT*

645 g

IMPACT LVL



WARRANTY

2 years

ENVIROHELL†



standard



With a focus on low-impact users, College Park combines some of its most popular and trusted technology in the new Terrain iW. This lightweight foot offers tri-axial movement on an Intelliweave composite base. Featuring anatomical ankle motion, flexible springs, and a durable heel wedge, the Terrain iW was designed for comfort and stability.

*26 cm with Enviroshell. ‡ The listing of codes with these products should not be construed as a guarantee for coverage or payment. Ultimate responsibility for the coding of services/products rests with the individual practitioner. †Caucasian, Tan, and Brown.

TERRAIN iW PART NUMBER

TN iW SIDE: L or R SHELL COLOR*: C, T, or B **EN** SIZE: 21-30 **D** FIRMNESS CATEGORY: 1-4

*(C)aucasian, (T)an or (B)rown. Example part number: TN iW LCEN 24 D3

TERRAIN iW FIRMNESS CATEGORIES

PATIENT WEIGHT	0-140lbs	141-180	181-220	221-300
SIZE	21-30cm		25-30cm	
LOW IMPACT	1	2	3	4

INCLUDED

Endo foot includes: sandal toe foot shell and CPI sock

ENVIROHELL PART NUMBER




Enviroshell included with every foot order






CS S SIDE: L or R SHELL COLOR**: C, T, or B SIZE: 21-30

** (C)aucasian, (T)an or (B)rown




HARDWARE AND TOOLS

	PART NUMBER	DESCRIPTION
	CP SSET	Super Set Screws (8pc) M8 x 14mm
	TFH 800	Foot Horn, Stainless Steel

	PART NUMBER	DESCRIPTION
	CP HK 4	Hex Key, 4mm

SHELLTREAD

	PART NUMBER (LEFT/RIGHT)	DESCRIPTION	ADDITIONAL INFORMATION
	CS TLN x CS TRN x	Shelltread (left or right), Neutral (indicate size)	Increases durability and prevents premature wear of foot shell.

SOCKS

	PART NUMBER	DESCRIPTION	ADDITIONAL INFORMATION
	CPI SS SM	CPI Sock, 21-25cm	Highly durable sock that acts as a noise barrier and protects the foot from debris.
	CPI SS LG	CPI Sock, 26-30cm	
	CPI SSK SM	CPI Sock, 6-pack, 21-25cm	
	CPI SSK LG	CPI Sock, 6-pack, 26-30cm	