















K3/K4 INTERNATIONAL PRODUCT DESK REFERENCE

T 586.294.7950 | F 586.294.0067 | college-park.com



	FOOT	MOUNTING	WEIGHT LIMIT	SIZES	IMPACT	BUILD HEIGHT	WEIGHT**	WTY	ENVIROHELL	HEEL HEIGHT
VERT SHOCK	 AVIATOR <small>FRESH WATER PROOF</small>	Endo (IP)	125 kg (21-25 cm) 150 kg (26-30 cm)	21-30 cm		16.4 cm	979 g	3		10 mm
	 VOYAGER <small>FRESH WATER PROOF</small>	Endo (IP)	125 kg (21-25 cm) 150 kg (26-30 cm)	21-30 cm		19.1 cm	1,174 g	3		10 mm
INTEGRATED SPRING TECHNOLOGY	 SOLEUS <small>WATERPROOF</small>	Endo (4-hole)	100 kg (21-24 cm) 113 kg (25-26 cm) 125 kg (27-30 cm)	21-30 cm		16.3 cm (21-24 cm) 18.3 cm (25-26 cm) 19.3 cm (27-30 cm)	667 g	3		10 mm
	 VELOCITY	Endo (IP)	113 kg (25-26 cm) 125 kg (27-30 cm)	25-30 cm		12.0 cm (25-26 cm) 13.2 cm (27-30 cm)	631 g	3		10 mm
TRI-AXIAL	 TERRAIN <small>WATERPROOF</small>	Endo (IP)	150 kg (21-25 cm) 166 kg (26-30 cm)	21-30 cm		8.6 cm (21-25 cm) 9.0 cm (26-30 cm)	815 g	3		10 mm
	 VENTURE	Endo (IP) Exo (ALX)	113 kg (21-24 cm) 125 kg (25-30 cm)*	21-30 cm		6.4-6.9 cm (IP) 12.7-13.2 cm (ALX)	585 g	3		10 mm
	 TRUSTEP	Endo (IP) Exo (ALX)	150 kg (Endo) 160 kg (Exo)	22-31 cm		8.7 cm (IP) 15.2 cm (ALX)	678 g	3		13 mm
	 TRIBUTE	Endo (IP)	100 kg	21-30 cm		6.1 cm	529 g	2		10 mm
CARBON+ TORSION	 ORBIT <small>FRESH WATER PROOF</small>	Endo (IP)	150 kg	21-30 cm		9.3 cm (21-25 cm) 9.8 cm (26-30 cm)	895 g	3		10 mm
CARBON	 HORIZON HD <small>FRESH WATER PROOF</small>	Endo (IP)	150 kg (21-25 cm) 166 kg (26-30 cm)	21-30 cm		5.8 cm	559 g	3		10 mm
K3 CURVED HYDRAULIC	 ODYSSEY_{K3}	Endo (IP)	125 kg	21-30 cm		7.5 cm (21-25 cm) 7.7 cm (26-30 cm)	749 g	3		10 mm
ADJUSTABLE HEEL	 ACCENT	Endo (IP) D.Pylon (DP)	100 kg (21-24 cm) 113 kg (25-28 cm)	21-28 cm		8.5-9.9 cm (IP) 19.8-20.2 cm (4-H) 20.3-20.6 cm (male) 21.1-21.5 cm (female)	608 g (IP) 761 g (DP) w/pylon	2		51 mm (range)
PEDIATRIC	 SCOUT <small>FRESH WATER PROOF</small>	Lam Adapter Pyramid Receiver Exo (ALX)	25 kg (13-15 cm) 45 kg (16-18 cm) 60 kg (19-22 cm)	13-22 cm		4.8-5.3 cm (Lam Adpt) 4.5-5.3 cm (Pyramid) 5.3-6.1 cm (Receiver) 11.0-11.8 cm (ALX)	183-360g [†]	1		6.1 mm (13-15cm) 6.4 mm (16-18cm) 9.7 mm (19-22cm)
	 TRUPER	Endo (IP) Exo (ALX)	45 kg (16-18 cm) 60 kg (19-21 cm)	16-21 cm		5.3-6.1 cm (IP) 11.0-11.8 cm (ALX)	215-314 g (16-21 cm)	1		6.3 mm

* Low impact only

**26 cm, including Enviroshell
†13-22 cm, including Pyramid Adapter and Enviroshell

Caucasian, Tan, Brown, Jet Black

	KNEE	K-LEVEL	WEIGHT LIMIT	WEIGHT	OVERALL HEIGHT	DOMES TO KNEE CENTER	FLEXION RANGE	WTY
MICRO-PROCESSOR	 ICON <small>NEW!</small> <small>IP68 FRESH WATER PROOF</small>	K2-K4	150 kg	1530 g	26.7 cm	1.7 cm	130°	3 (3 year extended warranty available)
HYDRAULIC	 CAPITAL <small>FRESH WATER PROOF</small>	K3-K4	150 kg	1010 g (pylon) 1020 g (pyramid)	23.0 cm	1.7 cm	130°	3
PNEUMATIC	 VICTOR	K3	125 kg	710 g	20.4 cm	2.9 cm	145°	3

K2 INTERNATIONAL PRODUCT DESK REFERENCE

T 800.728.7950 586.294.7950 | F 800.294.0067 586.294.0067 | college-park.com



	FOOT	MOUNTING	WEIGHT LIMIT	SIZES	IMPACT	BUILD HEIGHT	WEIGHT**	WTY	ENVIROSHIELD	HEEL HEIGHT	
K2 CURVED HYDRAULIC		ODYSSEY_{K2}	Endo (IP)	100 kg (21-24 cm) 136 kg (25-30 cm)	21-30 cm		7.8 cm (21-24 cm) 8.4 cm (25-30 cm)	744 g	2		10 mm
		TERRAIN iW	Endo (IP)	100 kg (21-24 cm) 136 kg (25-30 cm)	21-30 cm		6.4 cm (21-24 cm) 6.9 cm (25-26 cm) 7.1 cm (27-30 cm)	645 g	2		10 mm
LOW COST		CELSUS	Endo (IP)	100 kg (21-24 cm) 136 kg (25-30 cm)	21-30 cm		5.7 cm (21-24 cm) 6.4 cm (25-26 cm) 6.6 cm (27-30 cm)	526 g	3		10 mm
		TRES	Endo (IP)	136 kg (21-24 cm) 150 kg (25-30 cm)	21-30 cm		5.7 cm (21-24 cm) 6.4 cm (25-26 cm) 6.6 cm (27-30 cm)	500 g	2		10 mm
		BREEZE	Endo (IP)	100 kg (21-25 cm) 125 kg (26-30 cm)	21-30 cm		6.2 cm (21-25 cm) 6.5 cm (26-30 cm)	523 g (25 cm)	2		10 mm
K2 VACUUM		CELSUS iVT	Endo (IP)	100 kg (21-24 cm) 136 kg (25-30 cm)	21-30 cm		5.7 cm (21-24 cm) 6.4 cm (25-26 cm) 6.6 cm (27-30 cm)	652 g	2		10 mm
		TRES iVT	Endo (IP)	136 kg (21-24 cm) 150 kg (25-30 cm)	21-30 cm		5.7 cm (21-24 cm) 6.4 cm (25-26 cm) 6.6 cm (27-30 cm)	626 g	2		10 mm
STUBBIES		SIDEKICKS	Endo (IP)	125 kg	ONE SIZE		5.8 cm	422 g	2	No Enviroshell	0 mm

**26 cm, including Enviroshell

Caucasian, Tan, Brown

	KNEE	K-LEVEL	WEIGHT LIMIT	WEIGHT	OVERALL HEIGHT	DOMES TO KNEE CENTER	FLEXION RANGE	WTY	
POLYCENTRIC		AURA	K2	136 kg	665 g	16.1 cm	1.7 cm	135°	2
FRICTION BRAKE		GUARDIAN	K2	125 kg	533 g	13.1 cm	2.6 cm	145°	2

	LINER	TYPES	PROFILES	SIZES	OVERALL LENGTH	THICKNESSES	WTY	
TPE GEL		STERLING LINER	Cushion Locking	Uniform	S, M, M+, L, L+, XL	41.9 cm	3 mm, 6 mm, 9 mm	6 months
		STERLING SLEEVE	Suspension	Uniform	35, 41, 49, 55, 62	41.9 cm	5.5 mm at wave	90 days

ALSO AVAILABLE



ADVANCE YOUR EDUCATION

College Park Campus offers free education, training, and certification opportunities for prosthetists and other healthcare professionals.

Enroll today:
www.college-park.com/education



ENDO COMPONENTS

College Park components are highly engineered for noise abatement and strength.



TORSION ADAPTERS

Made of titanium and stainless steel, this adjustable component offers up to 20° of internal and external rotation.



BE THE FIRST TO KNOW

Sign up for the email newsletter to receive our latest news, events, promotions, and product releases!



LET'S CONNECT

We invite you to follow College Park on social media for helpful resources, inspiring stories, and more.

